ardo natal support

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Basel, 16.th of January 2012-01-16|

Dear Dr. Schlieper

As discussed before, we would like to inform you about our experience concerning birth gels used during deliveries.

We started to use birth gel, named Dianatal in 2008, while we were participating in a scientific study. This study showed that birth gels may facilitate the delivery, especially the second stage of labor. Since one year ago we collected experiences with the "Natalis" birth gel, which was developed by your company. We noticed that the Natalis birth gel has very similar physical characteristics as Dianatal. It can cause a comparable facilitation of the second stage of labor. In addition, Natalis offers further advantages compared with Dianatal. The production of Natalis is based on natural ingredients. Dianatal, however, contains polyacrylic acid as basis. Furthermore, the essentially softer application catheter of Natalis should be mentioned.

Detailed information of the Natalis birth concept was kindly provided by you. In our opinion the Natalis birth concept including Prenatal, Antenatal and Postnatal treatments offers a significant advantage to every woman in labour during natural birth. Essentially scientific data were considered in the birth concept of Natalis:

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Accredited by the European Board and College of Obstetrics and Gynaecology (E.B.C.O.G.)

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- Protection of the perineum and prophylaxis
- Facilitation of the delivery
- Reduction of the risk of infection for the mother after giving birth and in childbed

From our point of view applying the Natalis birth concept brings a positive impact into the field of natural birth.

Every pregnant woman who wants to give birth naturally and has no contraindication, may consider to use the birth concept of Natalis.

We are looking forward to continue our agreeable cooperation.

Yours sincerely

PD Dr. med. Olav Lapaire
University Women's Hospital Basel

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Obstetric gel shortens second stage of labor and prevents perineal trauma in nulliparous women: a randomized controlled trial on labor facilitation²

Andreas F. Schaub^a, Mario Litschgi^b, Irene Hoelsli^c, Wolfgang Holzgreve^c, Ulrich Bleul^d and Verena Geissbühler^e

^aHappy Child Birth Holding AG Basel, Switzerland

^bWomen's Hospital, Kantonsspital Schaffhausen, Switzerland

^cDepartment of Obstetrics, University Hospital Basel, Switzerland

^dClinic of Reproductive Medicine, Vetsuisse-Faculty, University Zurich, Switzerland

^eKantonsspital Frauenfeld, Switzerland

J. Perinat. Med. 36 (2008) 129-135

Objective: To determine whether the obstetric gel shortens the second stage of labor and exerts a protective effect on the perineum.

Method: A total of 251 nulliparous women with singleton low-risk pregnancies in vertex position at term were recruited. A total of 228 eligible women were randomly assigned to Group A, without obstetric gel use, or to Group B, obstetric gel use (Dianatal), i. e., intermittent application into the birth canal during vaginal examinations, starting at the early first stage of labor (prior to 4 cm dilation) and ending with delivery.

Results: A total of 183 cases were analyzed. For vaginal deliveries without interventions, such as C-section, vaginal operative procedure or Kristeller maneuver, obstetric gel use significantly shortened the second stage of labor by 26 min (30 %) (P = 0.026), and significantly reduced perineal tears (P = 0.024). First stage of labor and total labor duration were also shortened, but not significantly. Results did not show a significant change in secondary outcome parameters, such as intervention rates or maternal and newborn outcomes. No side effects were observed with obstetric gel use.

Conclusion: Systematic vaginal application of obstetric gel showed a significant reduction in the second stage of labor and a significant increase in perineal integrity. Future studies should further investigate the effect on intervention rates and maternal and neonatal outcome parameters.

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Use of obstetric gel in nulliparous pregnant women: Maternal and neonatal outcomes¹

Burçhan Aydiner^a, Hüseyin Kiyak^a, Fatih Mete^b, Ali Ekiz^a, Ibrahim Polaa^a, Ali Gedikbasi^a

^aGynecology and Obstetrics Clinic, Kanuni Sultan Süleyman Training and Research Hospital, Health Sciences University, Istanbul, Turkey

^bPediatrics Clinic, Kanuni Sultan Süleyman Training and Research Hospital, Health Sciences University, Istanbul, Turkey

Perinatal Journal 2017;25(3): 127-132

Objective: The study aims to investigate the effects of obstetric gel application on the first and second stages of labor in nulliparous pregnant women, the presence of its protective effect on perineum and fetal and maternal outcomes associated with this method.

Methods: The nulliparous pregnant women who admitted to our hospital between 37 and 41 weeks of gestation were randomized into two groups, which were those applied gel (Natalis*) and those not applied. During their labors, pregnant women were not administered any labor

induction or augmentation; they underwent routine amniotomy, but routine episiotomy was not performed during labor. Both groups were compared in terms of perineal laceration, episiotomy need, periods of labor stages, laboratory data, and newborn data.

Results: A total of 200 nulliparous pregnant women (Group I: n=102, spontaneous follow-up; Group II: n=98, gel-applied) who were complying with study criteria were evaluated. No difference was found in maternal demographic data, laboratory data and newborn results. In gel-applied patients, the first and second stages of labor were shorter and they were statistically significant. In terms of episiotomy procedure / need, there was no difference between the groups. In terms of cases where laceration did not occur and mild perineal lacerations (1st and 2nd degrees), the results were relatively better in gel-applied pregnant women, and no statistical difference was found in advanced perineal lacerations (3rd and 4th degrees).

Conclusion: Use of obstetric gel during labor in nulliparous women provides a significant reduction in labor stages and there are fewer lacerations.

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Effects of obstetric gel on neonatal and maternal outcomes in nulliparous women⁴

Dr Ghazala Amin^a, Dr Uzma Zahoor^b

^aSulaiman Al Habib Medical Group, Riyadh

^bDepartment of Gynecology, Takhussisi Branch, Sulaiman Al Habib Medical Group, Riyadh

Indo Am. J. P. Sci, 2019; 06(04)

Objective: The purpose of this study was to determine the effects of obstetric gel application on pregnant nulliparous women in the first and second stage of labor, to evaluate on the perineum with its protective effect, and the maternal and fetal outcomes related with this procedure.

Study Design: A prospective and randomized controlled study.

Place and duration: In the Obstetrics and Gynecology department of Sulaiman Al Habib Medical Group Hospital, Riyadh for one year duration from December 2017 to December 2018.

Methods: Into two groups; between 37 and 41 weeks of gestation nulliparous women admitted to our hospital were randomly divided. The pregnant women were not given any induction during their labour; Normal amniotomy was performed as routine, but during delivery routine episiotomy was not done. One Group with obstetric gel, the control group without. In both groups, episiotomy need, perineal laceration, laboratory data, newborn data and birth stages were compared.

Results: 200 total pregnant nulliparous females (n = 102 in Group I, in Group II: n = 98 spontaneous follow-up, treated with gel (Natalis*)) who met the criteria of the study were evaluated. There was no difference in laboratory data, neonatal results and maternal demographic data. The 1st and 2nd stages of labor were statistically significant and shorter. In terms of need for episiotomy/ procedure; there was no variation between the groups. As for the cases with mild perineal lacerations (first and second degree) and no fetal heart, the outcomes were better relatively in pregnant women with gel and there was no statistically significant variation in perineal lacerations at advanced stage (third and fourth grades).

Conclusion: In nulliparous women, during delivery the use of obstetric gels leads to a significant reduction in the stage of labor and less tearing.

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Obstetrical gel use to shorten labor and prevent lower genital tract trauma⁵

Elie Hobeika
American University of Beirut Medical Center

Published on ClinicalTrials.gov. 2019

Recent literature shows that birth injury is associated with postpartum pelvic floor dysfunction (pelvic organ prolapse and fecal and urinary incontinence). Prolonged labor, namely during the 2nd stage, is one of the main obstetric risk factors responsible for anal sphincter rupture and fecal incontinence. In addition, it is associated with increased maternal and neonatal morbidities including increased risk of lower genital tract lacerations. In an effort to shorten labor and decrease lower genital tract trauma many techniques have been investigated.

Objective: The objective of the study was to investigate whether the use of obstetric gel shortens the first and second stage of labor and exerts a protective effect on the lower genital tract. Neonatal and maternal morbidities was also assessed.

Study Design: The study design was a randomized controlled trial of 2 groups, where the patients presenting for vaginal delivery were randomly assigned to either:

- Group 1 who received the standard care during labor and delivery
- Group 2 who received the standard care during labor and delivery with the vaginal application of the obstetrical gel.

The goal of this randomized controlled study was to compare the length of the first and 2nd stage of labor and the lower genital tract integrity in the 2 groups of patients.

Methods: A specially-designed sterile obstetric gel (Natalis*) was used for this trial. The obstetric gel is a birth gel with no pharmacologic effects that has a purely physical activity. It contains propylene glycol, glycerol, sodium chloride, xanthan gum, hydroxyethylcellulose and water.

Starting with the first vaginal examination the obstetric gel was used. After each vaginal examination, 3-5 ml of obstetric gel was introduced into the vaginal birth canal in the area in front of the child using the sterile obstetric gel applicator without any manipulation or massaging. Additional obstetric gel was added 15-30 minutes after rupture of the membranes.

Primary Outcome Measures:

- Reduction in the length of the second stage of labor (Time Frame: Second stage of labor, a reduction to 39.9 minutes from a calculated average of 66.5 minutes)
- 40% reduction in the duration of the second stage of labor with a calculated average duration of 66.5 min. Thus a reduction to 39.9 minutes for every patient recruited.

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