## INSTRUCTIONS FOR USE

### Amniotic Membrane Perforator

en

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### Amnionmembran Perforator

de

## MODE D'EMPLOI

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it

#### Intended Use

The amniotic membrane perforator is an instrument used for the pregnant woman who needs artificial rupture of the amniotic membrane.

#### How to use

1. The use of the amniotic perforator must be performed by a qualified healthcare professional.
2. The perforator is sterile and must be opened in the delivery room.
3. Do not use if packaging is damaged.
4. The front hook is used to puncture; the other end is the handle. With the hook pointing upwards, hold the handle and insert the perforator into the cervix.
5. When the uterus contracts, press the hook against membrane, gently draw the instrument towards you, catching the membrane and rupturing it.
- Standard requirements of hospital operation must be used during the whole operation process.
6. After use, dispose of the perforator according to local medical waste requirement.

**Size:** L (265±1) mm, H (6.2±0.1) mm

**Sterilization method:** EO Sterilization

**Use by date:** Five years from the sterilization date.

**Material:** ABS (Acrylonitrile Butadiene Styrene)

#### Storage requirements:

Packaged products should be stored in an area at room temperature, normal humidity, non-toxic substances and well ventilated.

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Holding Corp. GmbH (Europe)
Hantech Medical Device Co., Ltd.
Ardo Medical AG